

# PRODUCT DESCRIPTION

Trolley for 3 litres lightweight cylinder LIV

Linde Gas  
Therapeutics

AGA



Complete trolley  
325 397 691



Bed hanger version  
325 197 566

This trolley has many combined functions. It enables many variations of applying equipment. The trolley exists in two models, complete trolley and bed hanger version.

The complete trolley has solid wheels and a strong robust frame. The trolley can either stand on its own or be hanged on a bed.

The bed hanger version can be carried by hand or placed on a bed gable.



## Article Numbers

Denomination	Article Number
Complete trolley	325 397 691
Bed hanger version	325 197 566

## Features

### Complete Trolley

- Cylinder diameter 118 mm.
- Basket for storage.
- Variable length of handle.
- Three sides to apply equipment.
- Six plate holders.
- Bed hanger.
- Holder to carry by hand.
- Cannula holder.

### Bed Hanger Version

- Special designed holder for carrying and placing on bed gables.

To equip the trolley, please contact your supplier.

## Recycling

The product shall be recycled according to local regulations.

**In order to guarantee a continuous technical development GCE reserve the right to make alterations in design and specification without prior notice.**

För att säkerställa kontinuerlig teknisk utveckling, förbehåller sig GCE rätten till ändringar utan föregående meddelande.

Manufactured by:  
GCE Norden AB  
Box 21044  
SE-200 21 Malmö  
Sweden

Product description  
700 325 738